Roche CARDIAC IQC cobas®

[REF]	[↓]	[SYSTEM]
04880668190	1 Roche CARDIAC IQC high test strip	cobas h 232
	1 Roche CARDIAC IQC low test strip	

English Intended use
The Roche CARDIAC IQC test serves as a performance check for the optical system of the **cobas h 232** instrument from Roche Diagnostics.

Summary
The Roche CARDIAC IQC test consists of two re-usable control strips to ensure correct functioning of the instrument.
The Roche CARDIAC IQC low test serves to check the lower measurement range.
The Roche CARDIAC IQC high test serves to check the higher measurement range.
Performance checks of the optical system should ideally be carried out on a daily basis.
The **cobas h 232** instrument should be checked using the Roche CARDIAC IQC test according to the applicable local or national regulations.

Reagents

One test contains:
Gold-labeled anti-HSA antibodies
IQC low: 85 ng
IQC high: 126 ng

Precautions and warnings

- For *in vitro* diagnostic use.
- Exercise the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in accordance with local guidelines.
- Safety data sheet available for professional user on request.
- Do not touch or wipe the signal line area of the test strip.
- Do not apply any sample material to the test strip.
- Do not expose the test strip to sunlight.
- Protect the test strip from dust and sunlight.
- Close the container immediately after removing the test strip to protect the second test strip from dust and moisture.
- Remove the test strip (low or high) from the **cobas h 232** instrument directly after the measurement is performed and place it quickly into the container to protect it from dust and moisture.
- For information on the optimum ambient temperature for the **cobas h 232** instrument, consult the operator’s manual for the instrument.

Lot code

Every kit contains a lot-specific code chip. The instrument display prompts the user to insert the chip. To ensure that the code chip and test strip lot match, compare the lot number in the display with the number on the code chip. The code chip provides the **cobas h 232** instrument with all the required lot-specific information. An error message is displayed if the wrong code chip is inserted for a test strip lot.

Test procedure

- Bring the Roche CARDIAC IQC test strips to room temperature before starting the measurement.
- Switch on the **cobas h 232** instrument.
- Remove a test strip, like Roche CARDIAC IQC low, from the container.
- Close the container immediately with the desiccant stopper.
- Choose "QC Test" from the main menu.
- When instructed in the display, insert the test strip with a smooth movement while holding it in an horizontal position. Be careful not to scratch the barcode on the underside of the strip. This can also happen if you push the strip over the small plastic positioning bar during insertion.
- Measurement starts as soon as the test reaches the measuring position.
- The result is displayed after a few seconds.
- Place the measured test strip quickly into the container and repeat the measurement with the second Roche CARDIAC IQC test strip from the container (start from step 3). For more details on the procedure, refer to the operator’s manual.

Storage and stability

Unopened, at 2-30 °C	up to the stated expiration date
After opening	6 months
	Write the opening date on the container label.

Materials provided

- [REF]** 04880668190, Roche CARDIAC IQC
- 1 Roche CARDIAC IQC high test strip
- 1 Roche CARDIAC IQC low test strip
- 1 code chip

Materials required (but not provided)

- cobas h 232** instrument

Display of results

On completion of measurement, the **cobas h 232** instrument shows the result in the display. The time required for the Roche CARDIAC IQC test to display a value is approximately 20 seconds. The **cobas h 232** instrument determines the reflectance value of the signal and displays "Pass" or "Fail".

"Pass" means that the optical system of the cobas h 232 instrument is functioning properly in the tested measurement range. If "Fail" appears, the value is outside the confidence interval.

Values outside the confidence interval:

- Repeat the instrument check using a new Roche CARDIAC IQC test strip of the same type (low or high) from a new, unopened container.
 - If the results obtained with the new test strip from the new container fall within the confidence interval, the first test strip is termed defective and must be discarded.
 - If the problem persists, contact your local Roche Diagnostics Technical Support team. For further information, please refer to the appropriate operator’s manual for the analyzer concerned, and the Method Sheets of all necessary components.
- A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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© 2019, Roche Diagnostics
COBAS, COBAS H, ELECSYS, IQC and ROCHE CARDIAC are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.roche.com