
Micral-Test

cobas[®]

REF	11544039171	∇ 30
REF	11544039172	∇ 30
REF	11544039005	∇ 30

English

Intended use

30 test strips for the immunological, semi-quantitative in vitro determination of urinary albumin or evaluation by visual reading. For professional use only.

Summary

Urinary albumin excretion of 20-200 mg/L¹ is called microalbuminuria. Microalbuminuria is an early indication of renal and cardiovascular diseases, which are both characterized by persistent albuminuria. Detection of microalbuminuria can aid diagnosis and treatment of incipient nephropathy in persons with diabetes and hypertension.² In addition, microalbuminuria is a predictor in the general population for cardiovascular outcomes independent of other risk factors such as hyperlipidemia, hypertension or diabetes.^{3,4}

Test principle

The albumin present in the urine specifically binds with a soluble antibody-gold conjugate present on a zone of the test strip. Excess conjugate is retained in a separation zone containing immobilized human albumin. This allows only the conjugate-albumin immunocomplex from the sample to reach the detection zone. After 1 minute, the intensity of the color produced (white to red) is directly proportional to the albumin content in the urine.

Reagents

Each test contains per 1 cm² reactive paper area the following: Monoclonal antibodies against human albumin (immunoglobulin G) labelled with colloidal gold: 6 µg, fixed albumin: 9.5 µg.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

Do not touch test zone or remove white covering foil from test strip. Micral-Test strips contain albumin of human origin.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods used assays approved by the FDA or cleared in compliance with the European Directive 98/79/EC, Annex II, List A. However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.^{5,6}

Note. In the following situations the detection of microalbuminuria does not always yield information with regard to hypertension-induced or diabetes-related impairment of renal function: acute illness, urinary tract infection, in urine testing positive for protein, nitrite, leukocytes or erythrocytes (e.g. with Combur¹⁰ Test strips), pregnancy, if physical exercise is done while urine is collecting in the bladder, if the metabolism is severely out of control or if albumin of post-renal origin is present.

If proteinuria is confirmed (e.g. a protein concentration of > 20 mg/dL or 200 mg/L or 0.2 g/L is measured with Combur¹⁰ Test strips), it is not usually necessary to screen for microalbuminuria. When comparing colors, always use the color scale on the container from which the test strips were taken.

The stopper of the test strip vial contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

Reagent handling

Test strips are ready for use.

Storage and stability

The test strips are stable up to the expiration date specified on the box, when stored in the original container.

Stability	
unopened at 2-8 °C	up to the stated expiration date

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Use only clean, well-rinsed vessels to collect urine.

Do not add urine preservatives.

Sample material: The first morning urine voided after rising is recommended as specimen material.² Physical activity can increase the excretion of albumin. If a random urine sample is used the albumin concentration is likely to be slightly higher than when using first morning urine.

⁷ Turbidity of the urine does not affect the test result.

Sample storage: If the urine is not tested within 3 days, it should be stored in a refrigerator (at 2-8 °C). Urine that has been refrigerated (for max. 2 weeks) must first be brought to at least 10 °C.

Diagnosis or therapy should never be based on one test result alone but should be established in the context of all other medical findings. In doubtful cases, it is therefore advisable to repeat the test after discontinuation of the medication.

Materials provided

For details see material table in header section.

Materials required (but not provided)

A vessel for collection of urine and a timepiece showing seconds.

Assay

For optimum performance of the assay follow the directions given in this document.

1. Take a test strip out of the vial. Close the vial again with the original desiccant stopper immediately after removal of the strip.

2. Place the test strip in the urine such that the fluid level is just between the 2 black bars (see arrows, Fig. 1), making sure that it does not touch the side of the vessel in the process. Withdraw the test strip after 5 seconds and place it across the top of the urine vessel.

3. After 1 minute compare the color of the test pad above the inscription "Micral" with the color scale on the test strip container label (Fig. 2). If color development is slightly uneven, it is the average color that counts. Comparison of the reaction color with the color scale is possible for another 5 minutes, as the color is stable for that period.

Evaluation: A wet test pad indicates that the reaction has come to an end. If the test pad is still dry after 1 minute despite correct immersion depth and duration, check the color development

after another 1 or 2 minutes. If the test pad is still dry, repeat the test with a new test strip paying careful attention to the correct depth and time of immersion.

Reaction colors lighter than the color block corresponding to approximately 20 mg/L albumin indicate a physiological urine albumin concentration (reference range). The screening result is positive when at least 2 of the 3 morning urines tested produce a reaction color corresponding to 20 mg/L albumin (threshold for microalbuminuria) or more. If the result is positive, note the concentration whose color block is closest to the test pad. In case it is unclear which color matches the test pad, choose a range, e.g. 20-50 mg/L or 50-100 mg/L. A positive screening result should be confirmed by nephrological examination. As albumin elimination is subject to physiological circadian variations⁷, tests should be performed on 2 different days, or on 3 different days if the results are contradictory.

Determination of albumin concentrations between 100 mg/L and 300 mg/L

In order to determine albumin concentrations between 100 mg/L and 300 mg/L, the urine sample can be diluted e.g. by mixing 1 part of urine with 5 parts of water. The determined albumin concentration should be assigned to the 50 mg/L color block. The original albumin concentration is then calculated by multiplying the result obtained by the dilution factor e.g. by 6.

Quality control

For quality control, use commercially available urine controls, or other suitable control material. The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits. Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

Therapeutic drugs and endogenous substances were tested for a potential interference to the Micral-Test. The Micral-Test was tested with negative urine samples and samples spiked to the first positive concentration range. The following therapeutic drugs and concentrations were tested:

Therapeutic drugs			
Substance	Maximum tested concentration	Substance	Maximum tested concentration
Acetaminophen	3000 mg/L	Levodopa	1250 mg/L
Ascorbic acid	4000 mg/L	Lisinopril	267 mg/L
Amoxicillin	6667 mg/L	Metformin	8500 mg/L
Biotin	1000 mg/L	Methyldopa	2000 mg/L
Cefoxitin	12000 mg/L	N-Acetylcysteine	200 mg/L
Furosemide	2000 mg/L	Ofloxacin	900 mg/L
Gabapentin	12000 mg/L	Phenazopyridine	300 mg/L
Gentamicin sulfate	400 mg/L	Salicylic acid	3000 mg/L
Ibuprofen	2500 mg/L	Tetracycline	500 mg/L

For therapeutic drugs testing within this range the following observations were made:

Therapeutic drugs	No interference up to	Effect above stated concentration
Ascorbic acid	400 mg/L	False-negative
Ofloxacin	100 mg/L	False-negative
Salicylic acid	100 mg/L	False-negative

Except of ascorbic acid, ofloxacin and salicylic acid, which lead to false-negative results, no interference by drugs has been found. Knowledge of the effects of drugs upon the test may not be complete. In doubtful cases, it is therefore advisable, as far as it is medically justifiable, to repeat the test after discontinuation of the medication. Endogenous substances were tested at abnormal high concentrations.

Endogenous substances			
Substance	Maximum tested concentration	Substance	Maximum tested concentration
Ammonium	25000 mg/L	Nitrite	110 mg/L
Calcium chloride	3000 mg/L	Bilirubin	1100 mg/L
Creatinine	15000 mg/L	Urea	200000 mg/L
α-D(+)-Glucose	50000 mg/L	Uric acid	1550 mg/L
Hemoglobin	750 mg/L	Urobilinogen	3000 mg/L
β-3-Hydroxybutyrate	4500 mg/L	pH	pH 5 and pH 9
Immunoglobulin G	5000 mg/L		

For endogenous substances testing within this range the following observations were made:

Endogenous substance	No interference up to	Effect above stated concentration
Ammonium	2500 mg/L	False-positive
Hemoglobin	300 mg/L	False-positive
β-3-Hydroxybutyrate	150 mg/L	False-positive
Uric acid	550 mg/L	False-positive
Urea	40000 mg/L	False-positive
Urobilinogen	120 mg/L	Elevated positive
Bilirubin	11 mg/L	No valid test result due to color interference

Potential cross-reactive urine proteins were tested for an interference to the Micral-Test. The Micral-Test was tested with negative urine samples and samples spiked to the first positive concentration range. The following substances and concentrations were tested:

Potential cross-reactive urine proteins			
Substance	Tested concentration	Substance	Tested concentration
Acidic α1-Glycoprotein	30 mg/L	α2-Macroglobulin	40 mg/L
α-Amylase (origin: pancreas)	3000 U/L	α1-Microglobulin	100 mg/L
α-Amylase (origin: salivary gland)	3000 U/L	β2-Microglobulin	300 mg/L
α1-Antitrypsin	100 mg/L	Retinol-binding protein	200 mg/L
Bence-Jones Protein	1800 mg/L	Tamm-Horsfall Protein	300 mg/L
Immunoglobulin A	400 mg/L	Transferrin	4 mg/L

The Bence-Jones protein leads to false-positive results and to elevated positive results. The Retinol-binding protein leads to elevated positive results.

Drug interferences are measured based on recommendations given in CLSI guidelines EP07 and EP37 and other published literature. Effects of concentrations exceeding these recommendations have not been characterized.

⚠ For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Sources of error

For reliable results on a concentration basis a normal fluid intake prior to testing is important (1.5 to 2 L fluids per day). A very low or very high fluid intake can lead to false-positive or false-negative results. False-negative results can be produced by residual quantities of strongly oxidizing cleaning agents in the urine collection vessel or if the test strip is immersed too far into the urine. Erroneous results are obtained if the strip is not immersed for the time specified, if the reaction color is not compared after the time specified or if the strip is allowed to touch the inside of the urine vessel. Urine samples which have been allowed to stand unrefrigerated for more than 3 days and show signs of bacterial decomposition (pH > 8) should no longer be used. If the urine specimen is colder than 10 °C the color reaction is diminished.

Expected and result values

Expected value < 20 mg/L⁸

Result values = neg., 20, 50, 100 mg/L

See color label on the test strip vial.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data are given below. Results obtained in individual laboratories may differ.

The values for neg. and pos. indicate the proportion of concordant negative or positive results. The values specified for the **limit of detection** (analytical sensitivity) are defined as the concentration of the analyte which leads to a positive result in > 90 % of the examined urines.

A **method comparison** was performed at 2-external sites using anonymized residual urine samples. The samples (n = 463 in total) were assayed by the Micral-Test strips and the Tina-quant ALB-T assay on a Roche Hitachi **cobas** c system as a quantitative comparison method (predicate). All concentration ranges were covered.

Parameter	Limit of detection	Method comparison ⁹⁾
Albumin	20 mg/L	neg.: 95 % (95 % confidence interval: 91.2 % - 97.6 %) pos.: 92 % (95 % confidence interval: 87.4 % - 94.6 %)

a) The values for neg. and pos. indicate the proportion of concordant negative or positive results.

Precision

Precision experiments comprised an assessment of repeatability (within-run precision) and intermediate precision. Repeatability was checked in 3 separate lots with 21 measurements each for the tested controls. In total there were 63 measurements performed per used control. The following results were obtained:

Repeatability			
Parameter	Control ⁹⁾	Result	Exact agreement
Albumin	Level 1	neg.	100 %
	Level 2	100 mg/L	100 %

b) Bio-Rad Liquichek Urinalysis Control

Intermediate precision was assessed over 20 days with 3 lots per day and four-fold measurements per used control. In total there were 80 measurements performed per used control and test strip lot. The following results were obtained:

Intermediate precision			
Parameter	Control ⁹⁾	Result	Exact agreement
Albumin	Level 1	neg.	100 %
	Level 2	50-100 mg/L	100 %

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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- References / Références bibliographiques / Referencias bibliográficas**
- McPherson RA, M.R.P., HENRY'S Clinical Diagnosis and Management by Laboratory Methods, 23rd edition. ISBN 9780323295680, 2017.
 - Hasslacher CH. Akt. Endokr. Stoffw. (1989); 10:60-63.
 - Amlöv J, Evans JC, Meigs JB, et al. Low-grade albuminuria and incidence of cardiovascular disease events in nonhypertensive and nondiabetic individuals: the Framingham Heart Study. Circulation (2005) Aug 16;112(7):969-75.
 - Klausen K, Borch-Johnsen K, Feldt-Rasmussen B, et al. Very low levels of microalbuminuria are associated with increased risk of coronary heart disease and death independently of renal function, hypertension, and diabetes. Circulation. 2004 Jul 6;110(1):32-5.
 - Occupational Safety and Health Standards: Bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
 - Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.
 - Microalbuminuria and Potential Confounders. Mogensen CE, Vestro E., Poulsen PL et al. Diabetes Care (1995); Vol.18(4).
 - Cembrowski G. Testing for Microalbuminuria: Promises and Pitfalls, Laboratory Medicine. 1990;21:491-496.

Symbols / Symboles / Símbolos
Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA; see https://usdiagnostics.roche.com for definition of symbols used):
/ Roche Diagnostics utilise les signes et les symboles suivants en plus de ceux de la norme ISO 15223-1 (pour les USA ; voir https://usdiagnostics.roche.com pour la définition des symboles utilisés) ; / Roche Diagnostics emplea los siguientes símbolos y signos adicionalmente a los indicados en la norma ISO 15223-1 (para los EE.UU.: consulte https://usdiagnostics.roche.com para la definición de los símbolos usados).

CONTENT	Contents of kit / Contenu du coffret / Contenido del estuche
SYSTEM	Analyzers/Instruments on which reagents can be used / Analyseurs/appareils compatibles avec les réactifs / Analizadores/instrumentos adecuados para los reactivos
REAGENT	Reagent / Réactif / Reactivo
CALIBRATOR	Calibrator / Calibrateur / Calibrador
➔	Volume after reconstitution or mixing / Volume après reconstitution ou homogénéisation / Volumen tras reconstitución o mezcla
GTIN	Global Trade Item Number / Code article international / Número mundial de artículo comercial

COBAS and MICRAL-TEST are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2019, Roche Diagnostics

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305

Mannheim

www.roche.com