Accutrend Control CH1 cobas®

REF	CONTENT	SYSTEM
11418289 190	1,5 mL	Accutrend GCT Accutrend Plus

English

Intended use

For performance checks of the Accutrend GCT or Accutrend Plus meter using

Accutrend Cholesterol test strips.

Suitable for self-testing.

Summary

Use Accutrend Control CH 1 solution to check the performance of Accutrend GCT or Accutrend Plus meter and the Accutrend Cholesterol test strips and their correct handling. Perform a performance check in the following situations:

- when learning how to measure cholesterol
- after changing batteries
- after opening a new pack of test strips and re-calibrating the meter
- after cleaning the meter
- if you have reasons to doubt the validity of a reading

How the test works

Each test strip has a test area containing detection reagents. When control solution is applied, a chemical reaction takes place and the test area changes colour. The meter records this change in colour and converts the measurement signal to the displayed result using the data previously entered via the code strip.

Important notes

Do not perform performance checks at temperatures below 18 °C or above 30 °C. Allow test strips and control solution just removed from the refrigerator to reach a temperature within this range before use. Use a fresh test strip if the drop of solution has not completely covered the test area, if several drops were applied or if a drop contained air bubbles. Do not touch the tip of the bottle.

Precautions and warnings:

For *in vitro* diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Safety data sheet available for professional user on request.

All components of the pack can be discarded along with household waste or, if used in a laboratory or doctor's office, according to relevant local guidelines.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

<div> </div> <div></div>	Warning
H317	May cause an allergic skin reaction.
Prevention:	
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves.
Response:	
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P333 + P313	If skin irritation or rash occurs: Get medical advice/attention.
P363	Wash contaminated clothing before reuse.
Disposal:	
P501	Dispose of contents/container to an approved waste disposal plant.
Product safety labeling primarily follows EU GHS guidance. Contact phone: all countries: +49-621-7590	

Storage and stability

Store at 2-25 °C up to the stated expiration date. Do not freeze. When stored at the conditions stated above, opened bottles of cholesterol control solution are stable for 3 months but not beyond the expiry date given on the pack. When opening the bottle for the first time, write the date onto the bottle label.

Sources of error

If a performance check with the Accutrend Control CH 1 solution produces a reading outside the stated range, please verify:

- that the cholesterol measurement was carried out correctly (see operator’s manual)
- that the test strip guide is clean (follow instructions for use)
- that the test strips and the control solution used were correctly stored (expiry date, storage conditions)
- that meter, test strips and control solution were stored at the correct temperature range of 18-30 °C during the measurement.

Preparing and performing the measurement

- Please read the information in this Method Sheet and in the operator’s manual for your meter.
- Carry out a performance check like you would do for a cholesterol measurement. Note: Use Accutrend Control CH 1 solution instead of blood.
- Remove one Accutrend Cholesterol test strip from the container and apply one drop of control solution to the test pad, moistening it completely.
- Hold the tip of the dropper bottle vertically downward and dispense one drop from a height of about 1 cm, without touching the test pad with the tip. Wipe the tip of the dropper bottle before and after use with a clean paper tissue or swab, and close the bottle tightly afterwards.
- Compare the result you get with the value sheet provided with the test strips. If the result is within the stated range, it shows that the meter and test strips are performing OK.
- If the control result is outside the stated range, check that you have followed the correct procedure (see section "Sources of error"). After eliminating all potential sources of error, repeat the performance check. If the control result is still outside the stated range after the second check, please contact your Roche Diagnostics representative.

Materials provided

- Accutrend Control CH 1: 1 x 1.5 mL of cholesterol control solution; sufficient for approximately 30 performance checks.

Materials required (but not provided)

- One of the above listed meters with operator’s manual
- Accutrend Cholesterol test strips
- Paper towel

Instructions for healthcare professionals

Applications: Perform quality control steps according to guidelines established for your practice or laboratory.

Reagents

Aqueous solution containing chemical process materials and extract of animal origin. The origin of the biological component is as follows: protein/lipoprotein serum (bovine). The concentration of the components is lot-specific and listed in the value sheet included in the kit.

Specific performance data

Reproducibility (day-to-day imprecision):

In a typical measurement series, the majority of the measured coefficients of variation (CV) was < 3.8 %.

Point sentence

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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